



www.gmmspa.com

# MAC D

Digital Mobile Radiography



Touch screen 8"  
console

Digital 19"  
workstation



## Post Processing

- Single image display.
- Multi-image overview.
- Dedicated LUT for any exam.
- Zoom and lens function
- Windowing /contrast adjustment.
- Level/brightness adjustment.
- Grey scale inversion.
- Image flips Left/right and Up/Down.
- Image rotation in 90° steps.
- Sharp and smooth spatial filters
- Electronic collimator with square/circular shutters.
- Text and markers annotation on the image.
- Graphic calculation of: angles and relative distances.
- DICOM – standard package

Different configurations are available within the MAC range with either Wi-Fi or wired detector : a wide series of advanced units tailored to the operator's actual needs.

## Generator Characteristics

|                            |                |
|----------------------------|----------------|
| Power kW                   | 32 kW          |
| Monoblock thermal capacity | 1000 kJ        |
| Focal spot mm              | 0.6/1.3 mm     |
| mAs range                  | 0.5 - 320      |
| kV range                   | 40 - 125       |
| Detector                   | Wi-Fi or wired |

## Applications

General radiology | Sports medicine | Emergency | Orthopaedics | Intensive care units | Operating rooms



## Distributor



**HORIZON MEDICA**  
MEDICAL EQUIPMENT SALES | MAINTENANCE | SERVICE

horizonmedica.com