



TRIUP TRF-100 Mobile X-Ray

INNOVATION ON THE MOVE.

X-Factor thinking has revolutionized portable imaging. The TRF-100 Mobile X-ray System lets you bring the key capabilities of a radiology room right to the patient at the bedside, in the OR, the ICU or EDU.



MOBILITY REDEFINED.



HORIZONMEDICA.COM Smart Services

A PARTNER
WHO CARES
AS MUCH AS
YOU DO.



A Community of Service and Support
Benet from all the advantages of our Customer
Success Network. We work continuously to
provide superb imaging performance, offer
you new innovations as your needs change,
and help you make the most of your budget
and resources. Our Customer Success
Network surrounds you with a dynamic team of
exceptional experts, with a Single Point of Entry for easy, customized access
to the right people in every situation.



You and your patients will benet from the expertise and best practices only
Horizonmedica.com can deliver, based on thousands of customer engagements
worldwide and our 100-year heritage in medical-imaging innovation.

Product Specifications:

Power Supply:

- Voltage AC : 220V ~ 22V
- Frequency : 50Hz ~ 0.5Hz
- Capacity : 8KVA
- Resistance : 1Ohm

Radiography:

- Voltage : 50-100kV
- Current : 15, 30, 60, 100mA
- Time : 0.08s-6.3s

Input Power :

- 5.92KVA

Performance :

- Distance between the focus of x-ray generator to the ground
UP and DOWN : 810-1790mm
- X-ray generator transverse retractility : 255mm
- X-ray generator LEFT and RIGHT : ~ 50 degrees
- X-ray generator move around level axis : ~ 180 degrees
- X-ray generator move around vertical axis : ~ 90 degrees

Dimension :

- 1442 (Length) x 690 (Width) x 1986 (Height) (mm)



Company:
Tel: +86-25-52334966
Fax: +86-25-52334165
Zip: 210019
Add: Floor #8, Building B, FengHuang Culture Plaza, JiangDong Middle Road#215 Nanjing China.



HORIZONMEDICA.COM
PT. HORIZON MULTI TEKINDO